

# Vacuum Experiment: Week 2

**AP 4018**  
**Columbia University**


# Week 2 Objectives

- Learn *more* about vacuum technology, instrumentation, and terminology
- Measure the pumping speed of the turbomolecular pump

# Introductory Reading

- Please see class homepage at: <http://sites.apam.columbia.edu/courses/ap4018y/> and
- Read background information about vacuum science and technology

*Very Good Introduction...*




[http://sites.apam.columbia.edu/courses/ap4018y/UHV\\_Seminar\\_Handbook.pdf](http://sites.apam.columbia.edu/courses/ap4018y/UHV_Seminar_Handbook.pdf)

# Overview of Week 2


- Adjust a controlled leak,  $Q_{\text{leak}}$
- Wait for the pressure to become steady, measure the pressure when the pumping rate exactly balances the leak rate.
- “Quickly” close the gate valve to the turbo pump, and allow the pressure to increase in the chamber. *This will be a linear increase.*
- Use the linear pressure increase to determine  $Q_{\text{leak}}$  and determine the pumping speed,  $S$ , of the turbo pump.
- Evaluate the pumping speed for different gases and equilibrium pressures.



Precision Leak Valve



Ionization Gauge (and Controller)



Read more online how ionization gauges work.

*Pressure indication:*

$$P_i = 10^{(V-11)} \text{ Torr or mbar}$$

where V is controller output in Volts



# Ionization Gauge Readings must be Corrected for Different Gas...

**Table 4-2** *Scale Factors for Use with Ionization Gauges.*

**(Do not use these scale factors for CONVECTRON Gauges)**

Gas	Scale Factor	Gas	Scale Factor
He	5.56	H <sub>2</sub> O	$8.93 \times 10^{-1}$
Ne	3.33	NO	$8.62 \times 10^{-1}$
D <sub>2</sub>	2.86	Ar	$7.75 \times 10^{-1}$
H <sub>2</sub>	2.17	CO <sub>2</sub>	$7.04 \times 10^{-1}$
N <sub>2</sub>	1.00	Kr	$5.15 \times 10^{-1}$
Air	1.00	SF	$4.00 \times 10^{-1}$
O <sub>2</sub>	$9.90 \times 10^{-1}$	Xe	$3.48 \times 10^{-1}$
CO	$9.52 \times 10^{-1}$	Hg	$2.75 \times 10^{-1}$

If the Stabil-Ion Gauge calibration is for a gas type other than N<sub>2</sub> (or air), we suggest placing a label near the first line of the display indicating the gas type or types used for calibration to prevent mix-ups.

(Why?)


# Procedure: Step 1

STEP #1

ADJUST LEAK VALVE TO A PARTICULAR EQUILIBRIUM PRESSURE.

$$V \frac{dP}{dt} \approx 0 \approx -SP + Q_{LEAK}$$

$$\therefore S \equiv \frac{Q_{LEAK}}{P_{EQUILIBRIUM}}$$




# Procedure: Step 2

STEP #2

CLOSE PUMP GATE VALVE. MEASURE LINEAR RATE OF CHANGE OF PRESSURE. (BE SURE TO USE EITHER IONIZATION GAUGE OR CONVECTION GAUGE BUT NOT BOTH FOR A SINGLE CONDITION)

$$V \frac{dP}{dt} = Q_{LEAK}$$

MEASURE  $Q_{LEAK}$ .  
CALCULATE S






# Procedure: Step 3

Repeat for different gas and different equilibrium pressure.



# Components of Your Vacuum System

- Oscilloscope can measure up to four signals simultaneously
- For Week #2, we record two:
  - Convector Analog Output
  - Ionization Gauge Analog Output





# Examples Available Today

Pumping Speed

